HI83314

Multiparameter Photometer with COD for Wastewater

with Digital pH Electrode Input

HI83314 benchtop photometer measures 10 different key wastewater quality parameters using 20 different methods that allow for multiple ranges and variations in chemistry for specific applications. The Chemical Oxygen Demand (COD) parameter is included for industrial and municipal wastewater treatment. The Phosphorous and Nitrogen parameters included are beneficial to municipal wastewater treatment customers that need to monitor their biological and chemical nutrient removal process. This photometer features an innovative optical system that uses LED's, narrow band interference filters, focusing lens and both a silicon photodetector for absorbance measurement and a reference detector to maintain a consistent light source ensures accurate and repeatable photometric readings every time.

To save valuable laboratory benchtop space, the HI83314 doubles as a professional pH meter with its digital pH/temperature electrode input. Now one meter can be used for both photometric and pH measurements.

Advanced optical system

 Innovative optical design that utilizes a reference detector and focusing lens to eliminate errors from changes in the light source and from imperfections in the glass cuvette.

• Backlit 128 x 64 Pixel Graphic LCD Display

- Backlit graphic display allows for easy viewing in low light conditions
- The 128 x 64 Pixel LCD allows for a simplified user interface with virtual keys and on-screen help to guide the user through use of the meter

Built-in Reaction Timer for Photometric Measurements

- The measurement is taken after the countdown timer expires.
- Countdown timer ensures that all readings are taken at the appropriate reaction intervals regardless of user for better consistency in measurements



Absorbance mode

- Hanna's exclusive CAL Check cuvettes for validation of light source and detector
- Allows for the user to plot concentration versus absorbance for a specific wavelength for use with user supplied chemistry or for teaching principles of photometry

• Units of Measure

 Appropriate unit of measure along with chemical form is displayed along with reading

Result Conversion

 Automatically convert readings to other chemical forms with the touch of a button

Cuvette Cover

 Aids in preventing stray light from affecting measurements

• Digital pH Electrode Input

- Measure pH and temperature with a single probe
- Good Laboratory Practice (GLP) to track calibration information including date, time, buffers used, offset and slope for traceability
- pH CAL Check alerts user to potential problems during the calibration process

 Space saving having a pH meter and photometer built into one meter

Data Logging

- Up to 1000 photometric and pH readings can be stored by simply pressing the dedicated LOG button. Logged readings are just as easily recalled by pressing the RCL button
- Sample ID and User ID information can be added to a logged reading using alphanumeric keypad

Connectivity

- Logged readings can be quickly and easily transferred to a flash drive using the USB-A host port or to a computer using the micro USB-B port
- Data is exported as a .CSV file for use with common spreadsheet programs

Rechargeable Battery

 Li-polymer rechargeable battery lasts for 500 measurements or 50 hours of pH measurement

Battery Status Indicator

· Indicates the amount of battery life left

Error Messages

- · Photometric error messages
- pH calibration messages include clean electrode, check buffer and check probe

Contact:

EASTERN ENERGY CO.,LTD. (HEAD OFFICE)
40/4 Vitoondumri Rd., Banbueng, Banbueng, Chonburi 20170
Tel: 0-3844-6117 Fax: 0-3844-6200 Email: info@ete.co.th



Specifications

Measurement Channels		5 x optical channels; 1 x digital electrode channel (pH measurement)	
Absorbance	Range	0.000 to 4.000 Abs	
	Resolution	0.001 Abs	
	Accuracy	±0.003 Abs (at 1.000 Abs)	
	Light Source	light-emitting diode	
	Bandpass Filter Bandwidth	8 nm	
	Bandpass Filter Wavelength Accuracy	±1.0 nm	
	Light Detector	silicon photocell	
	Cuvette Type	round, 24.6 mm diameter and 16 mm diameter	
	Number of Methods	128 max	
рН	Range	-2.00 to 16.00 pH (±1000 mV)*	
	Resolution	0.01 pH (0.1 mV)	
	Temperature Compensation	Automatic (-5.0 to 100.0°C; 23.0 to 212.0°F)*	
Temperature	Range	-20 to 120°C (-4.0 to 248.0 °F)	
	Resolution	0.1 °C (0.1 °F)	
Additional Specifications	pH electrode	digital pH electrode (not included)	
	Logging	1000 readings (mixed photometer and electrode); log on demand with user name and sample ID optional input	
	Display	128 x 64 pixel LCD with backlight	
	Connectivity	USB-A host for flash drive; micro-USB-B for power and computer connectivity	
	Battery Life	3.7 VDC Li-polymer rechargeable battery / >500 photometric measurements or 50 hours of continuous pH measurement	
	Power Supply	5 VDC USB 2.0 power adapter with USB-A to micro-USB-B cable (included)	
	Environment	0 to 50°C (32 to 122°F); 0 to 95% RH, non-condensing	
	Dimensions	206 x 177 x 97 mm (8.1 x 7.0 x 3.8 in.)	
	Weight	1.0 kg (2.2 lbs.)	

				LED (A nm) with Narrow Band	
Parameter Ra	ange	Resolution	Accuracy	Interference Filter	Method
Ammonia Low Range	0.00 to 3.00 mg/L (as NH ₃ -N)	0.01 mg/L	±0.04 mg/L ±4% of reading at 25 °C	@ 420 nm	Nessler
Ammonia Low Range (16 mm vial)	0.00 to 3.00 mg/L (as NH ₃ -N)	0.01 mg/L	± 0.10 mg/L or ± 5% of reading at 25 °C, whichever is greater	@ 420 nm	Nessler
Ammonia Medium Range	0.00 to 10.00 mg/L (as NH_3 - N)	0.01 mg/L	± 0.05 mg/L $\pm 5\%$ of reading at 25 °C	@ 420 nm	Nessler
Ammonia High Range	0.0 to 100.0 mg/L (a <mark>s NH</mark> ₃-N)	0.1 mg/L	± <mark>0.5 mg</mark> /L ±5% of reading at 25 °C	@ 420 nm	Nessler
Ammonia High Range (16 mm vial)	0.0 to 100.0 mg/L (as NH ₃ -N)	0.1 mg/L	± 1 mg/L or ± 5% of reading at 25 °C, whichever is greater	@ 420 nm	Nessler
Chlorine, Free	0.00 to 5.00 mg/L (as Cl ₂)	0.01 mg/L	± 0.03 mg/L $\pm 3\%$ of reading at 25 °C	@ 525 nm	DPD
Chlorine, Total	0.00 to 5.00 mg/L (as Cl ⁻)	0.01 mg/L	± 0.03 mg/L $\pm 3\%$ of reading at 25 °C	@ 525 nm	DPD
COD Low Range (16 mm vial)	0 to 150 mg/L (as O ₂)	1 mg/L	±5 mg/L or ±4% of reading @ 25 °C, whichever is greater	@ 420 nm	dichromate mercury-free
COD Medium Range (16 mm vial)	0 to 1500 mg/L (as 0 ₂)	1 mg/L	±15 mg/L or ±4% of reading @ 25 °C, whichever is greater	@ 610 nm	dichromate mercury-free
COD HR (16 mm vial)	0 to 15000 mg/L (as O_2)	1 mg/L	±150 mg/L or ±2% of reading @ 25 °C, whichever is greater	@ 610 nm	dichromate
Nitrate (16 mm vial)	0.0 to 30.0 mg/L (as NO ₃ -N)	0.1 mg/L	±1.0 mg/L or ±3% of reading at 25 °C, whichever is greater	@ 420 nm	chromotropic acid
Nitrite Ultra Low Range, Marine	0 to 200 μg/L (as NO ₂ - N)	1 μg/L	$\pm 10\mu g/L\pm 4\%$ of reading at 25 °C	@ 466 nm	diazotization
Nitrite Low Range	0 to 600 μg/L (as NO ₂ - N)	1 μg/L	$\pm 20\mu g/L\pm 4\%$ of reading at 25 °C	@ 466 nm	diazotization
Nitrite High Range	$0 \text{ to } 150 \text{ mg/L (as NO}_2^-\text{-N)}$	1 mg/L	± 4 mg/L $\pm 4\%$ of reading at 25 °C	@ 575 nm	ferrous sulfate
Nitrogen, Total Low Range (16 mm vial)	0.0 to 25.0 mg/L (as NO ₃ - N)	0.1 mg/L	±1.0 mg/L or ±5% of reading at 25 °C, whichever is greater	@ 420 nm	chromotropic acid
Nitrogen, Total High Range (16 mm vial)	0 to 150 mg/L (as N)	1 mg/L	±3 mg/L or ±4% of reading at 20 °C, whichever is greater	@ 420 nm	chromotropic acid
Phosphorus Reactive Low Range (16 mm vial)	0.00 to 1.60 mg/L (as P)	0.01 mg/L	±0.05 mg/L or ±4% of reading at 25 °C, whichever is greater	@ 610 nm	ascorbic acid
Phosphorus Reactive High Range (16 mm vial)	0.0 to 32.6 mg/L (as P)	0.1 mg/L	±0.5 mg/L or ±4% of reading at 25 °C, whichever is greater	@ 420 nm	vanadomolybdophosphoric acid
Phosphorus Acid Hydrolyzable (16 mm vial)	0 to 1.6 mg/L (ppm) (as P)	0.1 mg/L	±0.05 mg/L or ±5% of reading at 25 °C, whichever is greater	@ 610 nm	ascorbic acid
Phosphorus, Total Low Range (16 mm vial)	0.00 to 1.15 mg/L (as P)	0.01 mg/L	±0.05 mg/L or ±6% of reading at 25 °C, whichever is greater	@ 610 nm	ascorbic acid
Phosphorus, Total High Range (16 mm vial)	0.0 to 32.6 mg/L (as P)	0.1 mg/L	± 0.5 mg/L or $\pm 5\%$ of reading at 25 °C, whichever is greater	@ 420 nm	vanadomolybdophosphoric acid
Ordering Information HI83314-01 (115V) and HI83314-02 (230V) is supplied with sample cuvettes and caps (4 ea.), digestion vials (6), vial a cloth for wiping cuvettes, USB to micro USB cable connector, power adapter and instruction manual.					

Standards Contact:

EASTERN ENERGY CO.,LTD. (HEAD OFFICE)

W



HI83314-11 CAL Check Cuvette Kit for HI83399